## RILDResearch Institute for Liver Diseases (Shanghai) Co. Ltd

## Male Wistar Han (WH) Rat Hepatocytes

Product Number H		IP-WHDS-02M		Lot Number			SNDI	
Gender		Male		Num	Number of Animal in Pool		10	
Date of Preparation 2015.0			Storage Condition			-150° C		
(1) Cell Viability								
Specification							Results	
Post-thaw viability by trypen blue exclusion, >70%						90		
viable cell, >5 million/vial						5.4		
(2) Enzyme Characterization								
Phenotyping Reaction			Substrate Er		Enzyme	Enzyme activity (pmole/mg/min)		
			(μM) Mea		Mean		SD	
7-hydroxylation CoumarinGlucuronidation			30	0	116		17.2	
Testosterone 6β-hydroxylation			10	0	74.1		2.22	
(3) In vitro intrinsic clearance								
Probe Substrate	(μΜ)		T <sub>1/2</sub>	T <sub>1/2</sub> (min)		Clint(µL/min/10 <sup>6</sup> cells)		
Testosterone		1.0		4.0		173.89		
Verapamil		1.0		TBD		TBD		
120								
0 30 60 90 120 150 Time (min)								
Preparation by Junhua Fan		Fan	Signature		(	Innhua Fan		
Reviewed by Yi Sun		n	Signature			Junhua Fan Yi Sun Judhang Mus		
Approved by Yuchang		Mao	Signature			Judone Mus		

TBD: To be determined