

Beagle Dog Hepatocytes – DNEK

Product Number: HP-BGQ-01M

Lot No: DNEK

Storage Condition: -150°C

Test Results:

<u>Specification</u>	<u>Results</u>
> 70% post-thaw viability by trypan blue exclusion	71%
≥2.5 million viable cells	4.6

Lot Characterization Result:

Phenotyping Reaction	Substrate(μM)	Mean	SD
Acetaminophen	200	70.9	1.78
Testosterone 6β-hydroxylation	100	15.4	0.75
7-hydroxylation Coumarin Glucuronidation	300	149	4.49

Donor Information:

Species	Beagle Dog
Donor Number	1
Gender	Male
Age	9 month
Weight	11.5 kg

质量控制 _____ 李韶娟