

INVITROGRO™ HT Medium, serum-free, 225 mL
Product: S03317 | Lot: C19053B Expiration Date: 15 Dec 23

Storage: 4 °C

Test Results

	Specification	Results
Fill Volume	≥ 225 mL volume	225 mL
pH	7.40 – 7.60 at 20-30 °C	7.47

Add TORPEDO™ Antibiotic Mix at a ratio of 5.5 mL per 250 mL of medium. Medium expires 7 days after addition of antibiotic mix.

Results for this lot have been confirmed by Quality Assurance.

Caution: This product is for research use only. Do not use in animals or humans. This product has not been approved for any diagnostic or clinical procedures.