

Cynomolgus Monkey Liver S9 – RQRS

Cat. No.: LS-SXH-02F

Lot No.: RQRS

Unit Size: 0.5 mL

Storage Condition: -70°C

Protein Concentration: 20 mg/ml

Updated at: 2010-11

Lot Characterization Result

CYP450 Phenotyping Reaction	Enzyme Activity (pmol/ mg /min)		
	Substrate (μM)	Mean	SD
7-hydroxylation Coumarin Glucuronidation	300	429	21.9
Testosterone 6 β -hydroxylation*	200	617	59.7

Donor Information

Donor Number	7
Gender	Female
Naive	YES

质量控制

李韶娟