



**AROCLOR 1254-INDUCED  
MALE SPRAGUE-DAWLEY RAT LIVER S9  
PRODUCT NUMBER: M10002**

**Lot Number:** QEB

**Storage Conditions:** -70°C

**Test Results:**

Specification

20-26 mg/mL protein concentration  
total P450 concentration

Result

21.9 mg/mL  
0.481 nmol/mg

**Lot Characterization Results:**


Assay

ECOD: total rate of formation of 7-HC and metabolites

Result

3074 pmol/min/mg

**Number of donors in pool:** 20

  
\_\_\_\_\_

**Caution:** These products are for research use only. Do not use in animals or humans. These products have not been approved for any diagnostic or clinical procedures.