

Mouse Liver S9- OCRZ

Cat. No.: LS-XS-02F

Lot No.: OCRZ

Unit Size: 0.5mL

Storage Condition: -70°C

Protein Concentration: 20 mg/ml

Updated: 2011-02-25

Lot Characterization Result

Phenotyping Reaction	Enzyme Activity (pmol/mg/min)		
	Substrate (μ M)	Mean	SD
Testosterone 6 β -hydroxylation	40	383	2.34
7-hydroxylationcoumarin glucuronidation	300	318	13.2

Donor Information

Donor Number	30
Gender	Female
Age	8 周

质量控制

李韶娟