

Beagle Dog Liver S9- WOUF

Cat. No.: LS-BGQ-01F

Lot No.: WOUF

Unit Size: 0.5mL

Storage Condition: -70°C

Protein Concentration: 20 mg/ml

Updated:2011-02-25

Lot Characterization Result

Phenotyping Reaction	Enzyme Activity (pmol/mg/min)		
	Substrate (μ M)	Mean	SD
7-hydroxylation Coumarin Glucuronidation	300	758	40.2
Testosterone 6 β -hydroxylation	40	41.9	1.75

Donor Information

Species	Beagle
Donor Number	1
Gender	Female
Weight	10kg

质量控制

李韶娟