

Wistar Rat Liver S9 – KYUD
Cat. No.: LS-WDS-02M

Lot No. KYUD

Unit Size: 0.5 mL

Storage Condition: -70°C

Protein Concentration: 20 mg/mL

Updated at: 2012-2

Lot Characterization Result

Phenotyping Reaction	Enzyme Activity (pmole/mg protein/min)		
	Substrate (μM)	Mean	SD
Testosterone 6β-hydroxylation	40	343	42.3
7-hydroxylation Coumarin Glucuronidation	30	255	2.69

Donor Information

Species	Wistar rat
Donor Number	21
Gender	Male
Weight	200-250g

质量控制

李韶娟